



ACCESSORI COMPATIBILI PER APPARECCHIATURE COMPUMEDICS

- Model: 0510 Adult OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. 6-month warranty.
- Model: 0514PED Patented Pediatric size cTherm* Cannula Thermistor. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0580DPT2 Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0580LPT Lite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0529 High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.
- Model: 0510 Model: 0510PED Model: 0514 Model: 0514PED Model: 5512 Model: 0580 Model: 0580D Model: 0580L Model: 0528 Model: 0529 Model: 0570 Model: 0542 Model: 0540 Model: 0582-50 Model: 0582CO2-50 Model: 0582CO2-PED-50 Model: 0582PED-50 Model: 0582S-50 Model: 5582-40 Model: 0587-50 Model: 0589-50 Model: 0597-50 Model: 0599-25



Connector Cables

- 0574 18" (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.
- 0353 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) male phone connector.
- 0354 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) female phone connector.
- 0355 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) male phone connector.
- 0356 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) female phone connector.
- 0357 18" (45 cm) Dual 1 mm (keyhole) connector to RJ11 connector.
- 0358 18" (45 cm) Dual 1 mm (keyhole) connector to RJ10 connector.



REBEL

Healthcare Online


Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Accessori compatibili con apparecchiature Embla Embletta systems:

- Model: 8514cTherm Oronasal Airflow Thermistor for MediByte.
- Model: 0572 One-size-fits-most adult belt: adjustable from 24 inch to 48 inch (60 to 120 cm). One year warranty.
- Model: 0574 18 inch (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 8589-50 Ultima OroNasal Luer Lock Cannula (27 inch / 69 cm) - Qty 50.
- Model: 8597-50 PureFlow[®] OroNasal Cannula & Filter (27 inch / 69 cm) - Qty 50.

Accessori compatibili con apparecchiature Embla N7000 systems:

- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0572 One-size-fits-most adult belt: adjustable from 24 inch to 48 inch (60 to 120 cm). One year warranty.
- Model: 0573 18 inch (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.



Connector Cables

- 0574 18" (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.
- 0353 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) male phone connector.
- 0354 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) female phone connector.
- 0355 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) male phone connector.
- 0356 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) female phone connector.
- 0357 18" (45 cm) Dual 1 mm (keyhole) connector to RJ11 connector.
- 0358 18" (45 cm) Dual 1 mm (keyhole) connector to RJ10 connector.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Accessori compatibili con apparecchiature Grass systems:

- Model: 0510 Adult OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. 6-month warranty.
- Model: 0514PED Patented Pediatric size cTherm* Cannula Thermistor. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580 PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0580L PTLite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0529 High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 0530 5-Position Body Sensor - 6.5-foot cable length with male RCA connector.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.

Accessori compatibili con apparecchiature Nihon Khoden systems:

- Model: 0510 Adult OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. 6-month warranty.
- Model: 0514PED Patented Pediatric size cTherm* Cannula Thermistor. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580 PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0580L PTLite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0529 High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 053 2 5-Position Body Sensor - 6.5-foot cable length with 1.5 mm touchproof connector.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Accessori compatibili con apparecchiature Respironics/Philips

All options for Respironics Alice3 systems:

- Model: 0510A Alice 3 compatible Oronasal Thermistor. 7-month warranty.
- Model: 0580A Alice 3/4/5 compatible (includes 2 x 5582, 2 x 0589).
- Model: 0528A Alice 3 compatible with buckles (includes 0526-2, 0526-4).
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 053 25-Position Body Sensor - 6.5-foot cable length with 1.5 mm touchproof connector.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.

All options for Respironics Alice4 systems:

- Model: 0510A4 Alice 4 compatible Oronasal Thermistor. 7-month warranty.
- Model: 0580A Alice 3/4/5 compatible (includes 2 x 5582, 2 x 0589).
- Model: 0528A4 Alice 4 compatible with buckles (includes 0526-2, 0526-4).
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.

All options for Respironics Alice5 systems:

- Model: 0510A5 Alice 5 compatible Oronasal Thermistor. 7-month warranty.
- Model: 0575 18 inch (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0533 5-Position Body Sensor - 6.5-foot cable length for Alice 5. Variable Resistance.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Accessori compatibili con apparecchiature SANDMAN

All options for **Sandman SD20 systems:**

- Model: 0510SE Adult OroNasal Thermistor for Sandman SD20 amplifier. 7-month warranty.
- Model: 0514SE Patented cTherm* Cannula Thermistor for Sandman SD20 amplifier. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 0530M 5-Position Body Sensor - 6.5-foot cable length with 1/8 inch female phono connector.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.

All options for **Sandman SD32 systems:**

- Model: 0510 Adult OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. 6-month warranty.
- Model: 0514PED Patented Pediatric size cTherm* Cannula Thermistor. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580 PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0580L PTLite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0529 High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 0530M 5-Position Body Sensor - 6.5-foot cable length with 1/8 inch female phono connector.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Accessori compatibili con apparecchiature SensorMedics Cephalon systems

All options for **SensorMedics Cephalon systems**:

- Model: 0510 Adult OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. 6-month warranty.
- Model: 0514PED Patented Pediatric size cTherm* Cannula Thermistor. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580 PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0580L PTLite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0529 High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 0532 5-Position Body Sensor - 6.5-foot cable length with 1.5 mm touchproof connector.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.

Accessori compatibili con apparecchiature Stellate /Alpine

All options for **Stellate systems(Alpine)**

- Model: 0510 Adult OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. 6-month warranty.
- Model: 0514PED Patented Pediatric size cTherm* Cannula Thermistor. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580 PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0580L PTLite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0529 High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 0533 5-Position Body Sensor - 6.5-foot cable length for Alice 5. Variable Resistance.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.



REBEL


Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Accessori compatibili con apparecchiature Xltek systems

All options for **Xltek systems**:

- Model: 0510 Adult OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. 7-month warranty.
- Model: 0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. 6-month warranty.
- Model: 0514PED Patented Pediatric size cTherm* Cannula Thermistor. 6-month warranty.
- Model: 5512 Adult Oronasal Thermocouple with 1.5 mm safety pins. 6-month warranty.
- Model: 0580 PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- Model: 0580L PTLite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- Model: 0528 High Output with buckles & 1.5 mm safety pins (includes 0526-2, 0526-4). Alice 5 compatible.
- Model: 0529 High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4).
- Model: 0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts). 1-year warranty.
- Model: 0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck.
- Model: 0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face.
- Model: 0532 5-Position Body Sensor - 6.5-foot cable length with 1.5 mm touchproof connector.
- Model: 0582-50 Adult nasal cannula and filter - Qty 50.
- Model: 0582CO2-50 Adult CO2 Divided Cannula - Qty 50.
- Model: 0582CO2-PED-50 Pediatric CO2 Divided Cannula - Qty 50.
- Model: 0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.
- Model: 0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.
- Model: 0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.
- Model: 0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.
- Model: 0597-50 PureFlow OroNasal Cannula & Filter - Qty 50.
- Model: 0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.



Connector Cables

- 0574 18" (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.
- 0353 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) male phone connector.
- 0354 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) female phone connector.
- 0355 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) male phone connector.
- 0356 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) female phone connector.
- 0357 18" (45 cm) Dual 1 mm (keyhole) connector to RJ11 connector.
- 0358 18" (45 cm) Dual 1 mm (keyhole) connector to RJ10 connector.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Descrizione accessori



PureFlow® oronasal cannula & filter Model 0597

- Awarded best in class
- Best cannula for detecting both nasal and oral flow
- Combines nasal and oral flow into one channel - no internal turbulence
- Less position-sensitive than the other oral detecting cannulae
- Best filter on the market
- Single lumen, single-patient use, adult oronasal cannula & filter for any pressure sensor
- Also available for HST in a shorter version

PureFlow® Duo oronasal cannula & filters Model 0588



- Dual lumen, single-patient use, adult oronasal cannula & filters for PT2 Dual Pressure Sensor
- Best dual lumen cannula for detecting both nasal and oral flow on separate channels
- Maintains separate nasal and oral flow into two channels - no internal turbulence
- Less position-sensitive than the other oral detecting cannulae
- Best filters on the market



Nasal & oral luer lock cannula & safety filter Model 0589

- Single lumen, single-patient use, adult oronasal cannula & filter for any pressure sensor
- Combines nasal and oral flow into one channel - no internal turbulence



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com



Large bore nasal & oral cannula & filter Model 0587

- Single lumen, single-patient use, adult oronasal cannula & filter for any pressure sensor
- Combines nasal and oral flow into one channel - no internal turbulence



Adjustable nasal & oral luer lock cannula & filter Model 0599

- Single lumen, single-patient use, adult oronasal cannula & filter for any pressure sensor
- Combines nasal and oral flow into one channel - no internal turbulence



ValueFlow nasal cannula with male luer lock Model 5582

- Adult micro nasal cannula
- Available with or without filter
- Works with any pressure sensor



Large bore nasal cannula & safety filter Model 0582

- Adult nasal cannula
- Works with any pressure sensor



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com



Model 0582MCO₂

- Simultaneous CO₂ and pressure monitoring
- Male luer

Model 0582CO₂

- Simultaneous CO₂ and pressure monitoring
- Female luer

Which configuration do you need?

Cannulae and Accessories

Part # Description

0303 Male-to-Male Luer Lock Adapter - Qty 50.

0310 Pneumo Flow Sensor Kit.

0410 3.6 volt ½ AA lithium battery - Qty 3.

0582-50 Adult CO₂ Divided Cannula - Qty 50.

0582CO₂-50 Adult CO₂ Divided Cannula - Qty 50.

0582CO₂-PED-50 Pediatric CO₂ Divided Cannula - Qty 50.

0582PED-50 Pediatric Nasal Cannula & Safety Filter - Qty 50.

0582S-50 Adult Micro Nasal Luer Lock Cannula & Safety Filter - Qty 50.

0582S-3-50 Adult Nasal Cannula, Single Use, 3' (90 cm), male luer, no filter, 50 per case.

0582S-7-50 Adult Nasal Cannula, Single Use, 7' (210 cm), male luer, no filter, 50 per case.

0587-50 Large Bore Nasal & Oral Cannula & Safety Filter - Qty 50.

0589-50 Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 50.

0597-50 PureFlow OroNasal Cannula & Filter - Qty 50 - Best OroNasal Cannula on the market.

0588-50 PureFlow Duo OroNasal Cannula(tm) & Filters - Qty 50 - Dual Lumen - use only with PT2 Dual Pressure Sensor

0599-25 Adjustable Nasal & Oral Luer Lock Cannula & Safety Filter - Qty 25.

0583-50 Safety Filter - Qty 50.

0583Pes-50 Safety Filter (includes 0303 male-to-male luer lock adapters) - Qty 50.

0591-5 Esophageal Catheter - Qty 5.

5582-40 ValueFlow Nasal Cannula with male luer lock and filter - Qty 40.

5589-22 ValueFlow Plus Nasal + Oral Cannula with male luer lock and filter - Qty 22.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Airflow Sensor Technology



Airflow temperature sensor
Models 0510, 0510SE, 0510Ped

- Most comfortable design available
- No prongs to go up the nose
- Connects directly to the headbox
- Works with all PSG systems
- Longest warranty in the industry
- Warranty: 7 months



Patented cTherm Cannula Thermistor
Airflow temperature sensor
Models 0514, 0514A3, 05104A5, 0514Ped, 0514A5Ped

- Works well with any cannula
- Excellent thermal airflow and pressure signals with all cannulae
- No prongs to go up the nose
- Connects directly to the headbox
- Designed to comply with new guidelines [Learn more](#)



Airflow sensor
Model 5512

- Thermocouple technology
- Measures both oral and nasal airflow
- Two 1.5mm safety connectors
- Cable length 84 inches (210 cm)
- Warranty: 6 months



Airflow temperature sensor
Model B1221

- Ideal for use during CPAP titration
- No prongs to go up the nose
- Connects directly to the headbox
- Longest warranty in the industry
- Warranty: 6 months



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Which configuration do you need?

Thermistors

Part # Description

0510 Adult OroNasal Thermistor with 1.5 mm safety pins. Compatible with most PSG systems. 7-month warranty.

0510PED Pediatric (ages 1-10 yrs) OroNasal Thermistor with 1.5 mm safety pins. Compatible with most PSG systems. 7-month warranty.

0510A Alice 3 compatible Oronasal Thermistor. 7-month warranty.

0510A4 Alice 4 compatible Oronasal Thermistor. 7-month warranty.

0510A5 Alice 5 compatible Oronasal Thermistor. 7-month warranty.

0510BS Reusable BreathSensor cable compatible with NPB disposable airflow sensor. 7-month warranty.

0510SE Adult OroNasal Thermistor for Sandman SD20 amplifier. 7-month warranty.

0514 Patented cTherm* Cannula Thermistor with 1.5 mm safety pins. Compatible with most PSG systems. 6-month warranty.

0514A3 Patented cTherm* Cannula Thermistor for Alice 3 amplifier. 6-month warranty.

0514A5 Patented cTherm* Cannula Thermistor for Alice 5 amplifier. 6-month warranty.

0514PED Patented Pediatric size cTherm* Cannula Thermistor. Compatible with most PSG systems. 6-month warranty.

0514A5PED Patented Pediatric size cTherm* Cannula Thermistor for Alice 5 amplifier. 6-month warranty.

0514SE Patented cTherm* Cannula Thermistor for Sandman SD20 amplifier. 6-month warranty.

0515 Oral only Thermistor. Compatible with most PSG systems. 7-month warranty.

Thermocouple

Part # Description

5512 Adult Nasal/Oral Thermocouple with 1.5 mm safety pins. Compatible with most PSG systems. 6-month warranty.

B1221 Reusable oral only thermocouple ideal for use during CPAP titration. Compatible with most PSG systems. 6-month warranty.

*Alice is a trademark of Phillips.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Snore Sensor Technology



Snore Microphone Model 0540

- Electret condenser microphone
- Designed to detect higher frequencies of snoring
- Connects to headbox and senses snoring sounds
- Smaller than a dime, allowing comfortable placement
- Warranty: 6 months



Snore Sensor Model 0542

- Piezo technology
- Small and light
- Two 1.5 mm safety connectors
- Cable length 84 inches (210 cm)
- Warranty: 9 months



Q-Snor Kit Model 0546

- Quantitative snoring sound volume meter
- Comes with a Q-Snor microphone (Model 0545)
- Accurately and reliably record decibel sound level intensities
- True accurate pre and post-therapy study comparisons
- Two 1.5mm safety connectors
- Cable length 84 inches (210 cm)
- Warranty: 1 year (interface)



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com



Snore Microphone for MediByte® Model 8541

- Electret condenser microphone
- Compact and lightweight
- Locking touchproof connector
- Cable length 36 inches (90 cm)
- Warranty: 6 Months

Which configuration do you need?

Snore Sensors

Part # Description

0540 Snore Microphone - 6.5-foot cable length with 1.5 mm safety pins - place on face. Compatible with Most PSG Systems.

0542 Snore Sensor - 6.5-foot cable length with 1.5 mm safety pins - place on neck. Compatible with Most PSG Systems.

8541 Snore Microphone for MediByte® - 36 inch (90 cm) cable length with locking touchproof connector - Place on face.

Quantitative Snore System

Part # Description

0544 Q-Snor™ Interface.

0545 Q-Snor™ Microphone.

0546 Q-Snor™ Kit (includes Model 0544 Interface and Model 0545 Microphone).

Connector Cables

0574 18" (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.

0353 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) male phone connector.

0354 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) female phone connector.

0355 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) male phone connector.

0356 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) female phone connector.

0357 18" (45 cm) Dual 1 mm (keyhole) connector to RJ11 connector.

0358 18" (45 cm) Dual 1 mm (keyhole) connector to RJ10 connector.



REBEL


Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Effortless Effort Measurement

RIP Technology

New guidelines require the use of inductive plethysmography for effort measurement and BRAEBON's new Q-RIP™ effort system is the better solution! Our one-size-fits-most belt uses the correct X-Y axis for inductance measurement. The interface unit always includes the sum output and the unit is smart - it turns itself off after 14 hours. This premium high-quality product is durable and backed by a one-year warranty on belts and interface.



Q-RIP Respiratory Inductive Plethysmography Kit
Model 0570

- Superior quality RIP system
- Senses changes in respiratory effort
- Features 14 hour auto-off
- Detect chest and abdomen effort as well as SUM
- Kit includes the belts and all required cables
- Replaceable AA rechargeable or disposable batteries
- Warranty: 1 year



Q-RIP Belts
Models 0572, 0572L, 0572PED

- Premium Respiratory Inductive Plethysmography belts
- 0572 - One-size-fits-most adult belt (24" - 48")
- 0572L - Extra adjustable RIP belt (48" - 72")
- 0572PED - Pediatric adjustable RIP belt (12" - 24")
- Warranty: 1 year



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Piezo Technology

Our Piezo Respiratory Effort Belts are compatible with most polygraphs and computerized PSG systems and install quickly and easily. Connect them directly to your AC headbox - no batteries required!



Respiratory Effort Sensor Model 0528

- High output up to five times stronger than others
- Designed to reduce false paradoxical effort
- Connects directly to the head box
- Available for most PSG systems in pediatric and adult sizes
- Warranty: 1 year

Piezo belt models

- 0520: 1.5 mm safety pins with loops - 6.5' cable (includes 0523-3, 0523-5)
- 0522: 1.5 mm safety pins with buckles - 6.5' cable (includes 0526-3, 0526-5)
- 0528: High Output with buckles & 1.5 mm safety pins (Alice* 5 Compatible) (shown)
- 0528A: Alice* 3 compatible with buckles
- 0528A4: Alice* 4 compatible with buckles
- 0529: High Output with loops & 1.5 mm safety pins (includes 0523-2, 0523-4)



Piezo Accessories

- 0523-2: 2-foot Velstretch Loop Belt (orange)
- 0523-3: 3-foot Velstretch Loop Belt (black)
- 0523-4: 4-foot Velstretch Loop Belt (green)
- 0523-5: 5-foot Velstretch Loop Belt (red)
- 0526-2: 2-foot Velstretch Buckle Belt (orange)
- 0526-3: 3-foot Velstretch Buckle Belt (black)
- 0526-4: 4-foot Velstretch Buckle Belt (green)
- 0526-5: 5-foot Velstretch Buckle Belt (red)



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Which configuration do you need?

Q-RIP

Part # Description

0570 Q-RIP Starter Kit (Includes one 0571 interface plus two 0572 belts).

0571 Q-RIP interface with SUM Output. One year warranty. Uses replaceable AA rechargeable or disposable batteries.

0572 One-size-fits-most adult belt: adjustable from 24" to 48" (60 to 120 cm). One year warranty.

0572Ped Pediatric Adjustable RIP Belt: adjustable from 12" to 24" (30 to 60 cm). One year warranty.

0572L Extra Large Adjustable RIP Belt: adjustable from 48" to 72" (120 to 180 cm). One year warranty.

Connector Cables

0574 18" (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.

0353 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) male phone connector.

0354 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) female phone connector.

0355 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) male phone connector.

0356 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) female phone connector.

0357 18" (45 cm) Dual 1 mm (keyhole) connector to RJ11 connector.

0358 18" (45 cm) Dual 1 mm (keyhole) connector to RJ10 connector.

Piezo Respiratory Effort Sensors

Part # Description

0520 1.5 mm safety pins with loops - 6.5' cable. Compatible with most PSG Systems (includes 0523-3, 0523-5)

0522 1.5 mm safety pins with buckles - 6.5' cable. Compatible with most PSG Systems (includes 0526-3, 0526-5).

0528 High Output with buckles & 1.5 mm safety pins - Compatible with most PSG Systems including Alice 5 (includes 0526-2, 0526-4).

0528A Alice 3 compatible with buckles (includes 0526-2, 0526-4).

0528A4 Alice 4 compatible with buckles (includes 0526-2, 0526-4).

0528M High Output with buckles & 1.5 mm safety pins that connects with BRAEBON's PSG System.

0529 High Output with loops & 1.5 mm safety pins. Compatible with most PSG Systems (includes 0523-2, 0523-4).

5527 Respiratory Effort Sensor with 1.5 mm safety pins. Compatible with most PSG Systems.

Accessories

0523-2 2-foot Velstretch Loop Belt

0523-3 3-foot Velstretch Loop Belt

0523-4 4-foot Velstretch Loop Belt

0523-5 5-foot Velstretch Loop Belt

0526-2 2-foot Velstretch Buckle Belt

0526-3 3-foot Velstretch Buckle Belt

0526-4 4-foot Velstretch Buckle Belt

0526-5 5-foot Velstretch Buckle Belt

*Alice is a trademark of Phillips.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Differential Pressure Transducers

These compact battery operated pressure transducers use state-of-the-art-miniaturized technology to detect and reproduce waveforms associated with respiratory airflow and snoring. Developed by an experienced team of industry leaders, our transducers ensure superior quality and provide all you need to capture high quality and reliable recordings.



PTLite
Low gain differential pressure transducer starter kit
Model 0580L

- Includes 2 x 0589 and 3 x 0582S cannulae
- Ultima Pressure Sensor Interface - Model 0586L included
- Ultima Pressure Sensor - Model 0581L
- Up to three outputs: raw airflow, filtered airflow and snoring
- Warranty: 1 year



PT1
Low gain/High gain differential pressure transducer starter kit
Model 0580

- Includes 2 x 0589 and 3 x 0582S cannulae
- Ultima Pressure Sensor - Model 0581
- Compatible with most PSG systems - various connectors available
- Up to three outputs: raw airflow, filtered airflow and snoring
- Warranty: 1 year



PT2 Dual
Dual airflow differential pressure transducer
Model 0580D

- Includes 2 x 0582S and 2 x 0589
- Compatible with most PSG systems - various connectors available
- Up to six outputs: raw airflow, filtered airflow and snoring
- Warranty: 1 year



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com



Pes Kit Esophageal pressure sensor Model 0580Pes

- Includes 2 x 0582S, 2 x 0589, 2 x 0591, 2 x 0303
- Ideal for measuring endoesophageal pressure when using recommended catheter
- Warranty: 1 year



CPAP Titration Kit For titrating FPH CPAPs Model 0550

- Detect CPAP flow, leak, pressure and snoring
- Interfaces to your PSG system
- User-replaceable AA batteries
- Warranty: 1 year



Esophageal Catheter Model 0591-5

- Metal stylet
- Centimeter markings
- 5 Fr size
- Adult Single Patient Use
- Shipped sterile
- Sold in packs of five



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Piezo Pressure Sensor Model 5581



- Piezo technology
- Both Airflow and Snoring output signals
- Two pairs of 1.5 mm safety connectors
- Plugs directly into the headbox
- Warranty: 6 Months

Which configuration do you need?

Pressure Transducers

Part # Description

- 0580L PTLite Pressure Transducer AC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- 0580 PT1 Pressure Transducer AC/DC Starter Kit - (includes 2 x 0589 & 3 x 0582S cannulas).
- 0580A Alice 3/4/5 compatible (includes 2 x 5582, 2 x 0589).
- 0580D PT2Dual Airflow Pressure Transducer (includes 2 x 0582S, 2 x 0589).
- 0580Pes Esophageal Pressure Sensor (includes 2 x 5582, 2 x 0589, 2 x 0591, 2 x 0303). Requires DC input.
- 0591-5 Esophageal Catheter - Qty 5.

Titration Sensors

Part # Description

- 0550 CPAP Titration Starter Kit (incl. Model 0551, 0552, 0553, choice of four 18" adapter cables).
- 0551 CPAP Titration Pressure Sensor.
- 0552 Pressure Tubing Kit connects 0551 Titration Sensor to CPAP mask connector (Qty 5).
- 0553 Pitot Tubing Kit connects 0551 Titration Sensor to FPH Pitot Tube.

Connector Cables

- 0574 18" (45 cm) Dual 1 mm (keyhole) connector to two 1.5 mm touchproof connectors.
- 0353 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) male phone connector.
- 0354 18" (45 cm) Dual 1 mm (keyhole) connector to 1/8" (3.5 mm) female phone connector.
- 0355 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) male phone connector.
- 0356 18" (45 cm) Dual 1 mm (keyhole) connector to 3/32" (2.5 mm) female phone connector.
- 0357 18" (45 cm) Dual 1 mm (keyhole) connector to RJ11 connector.
- 0358 18" (45 cm) Dual 1 mm (keyhole) connector to RJ10 connector.

Pressure Sensor

Part # Description

- 5581 Piezo Pressure Sensor - conveniently detects both airflow and snoring without batteries.

*Alice is a trademark of Phillips.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Precise Body Position

Braebon DC sensors are smaller, smarter and easier to use than ever before. Accelerometer-based algorithms ensure signal accuracy.



Body Position Sensor

Models 0530, 0530M, 0531, 0532, 0533

- Accelerometer-based 5-position sensor
- Easily attaches to effort belt
- Two-year battery life
- Available for most PSG systems
- Warranty: 1 year

Signal Output

Body Position Voltage ± 0.1

Right Side 0.50 volts

Left Side 1.00 volts

Supine 0.25 volts

Prone 0.75 volts

Upright 0.00 volts

Note: 0533 uses variable resistance;
voltage values do not apply.

Which configuration do you need?

Body Position Sensors

Part # Description

0530 5-Position Body Sensor - 6.5-foot cable length with male RCA connector.

0530M 5-Position Body Sensor - 6.5-foot cable length with 1/8" female phono connector.

0531 5-Position Body Sensor - 6.5-foot cable length for MediPalm®.

0532 5-Position Body Sensor - 6.5-foot cable length with 1.5 mm touchproof connector.

0533 5-Position Body Sensor - 6.5-foot cable length for Alice 5.

*Alice is a trademark of Phillips.



REBEL


Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com



The perfect pair

Experience the enhanced comfort and performance of the innovative PureFlow[®] cannula. For optimum results, pair it with our patented cTherm[®] cannula thermistor to record airflow and nasal pressure simultaneously.



Separate chambers

Get the best possible results with an innovative cannula featuring separate chambers for oral and nasal flow. This separation increases the accuracy of true hypopnea detection, improving the results for interpretation.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com



Large oral capture

Increase oral flow to obtain incredible signal quality with the largest oral capture area on the market. This translates to a greater tolerance to initial placement, greater tolerance to patient movement and higher confidence in success.



Larger filter, ready to wear

A better, larger filter minimizes resistance to flow, producing higher quality data capture and delivers a more comfortable experience.

To make study preparation easier, the PureFlow[®] does not require any trimming or adjusting of the oral tube.



REBEL

Healthcare Online

Rebel snc
Via Koch, 31/A – 43123 Pilastrello – Parma (PR)
Tel. 0521 346603 - Fax +39 0521 641804
info@rebelmed.com – www.rebelmed.com

Disposable Airflow



The newest couple

Taking comfort and performance to a new level, the innovative 0518 Disposable Airflow Sensor pairs with the PureFlow® cannula to deliver simultaneous capture of airflow and nasal pressure.

The 0518 Disposable Airflow Sensor requires the BRAEBON connector cable model 0510BS.



Thin and powerful

Maintain the same incredible signal quality as our reusable cTherm® cannula thermistor in a super thin, disposable form

Nasal Cannula. No problem.

The same form-fitting 0518 easily fits a nasal-only cannula by trimming the extra side wraps.
(No image)

Great on its own

Using incredibly comfortable adhesive, the standalone 0517 Disposable Airflow Sensor offers the same comfort and performance as the 0518 without the use of a separate cannula.
(No Image)